



INSTRUCTIONS FOR USE

Protective Disposable Respirator With Valve

MODEL: **ALIVE PPE 9020**

PURPOSE: Alive PPE protective disposable respirator with valve is designed to release hot, humid exhaled breath quickly, helping to reduce heat build-up and moisture inside the face piece. Protective disposable respirator with valve made with;

-outer layer is spunbond polypropylene 60 g/m² (+ - 2 gr)

-middle layer is PA6 2 g/m² (+ - 0,5 gr)

-inner layer is spunbond polypropylene 60 g/m² (+ - 2 gr)

PROPERTIES:

Nose Bridge Wire: PE plastic covered copper wire 1 mm width, 95 mm length

Elastic band: PES with EA band (Ear band) 4 mm width, 210 mm length/ear

Exhalation Valve: Designed to release hot, humid exhaled breath quickly, helping to reduce heat build-up and moisture inside the facepiece.

LIMITATIONS:

Exposure to certain chemicals or high concentrations may require higher protective properties through the material or construction properties of the respirator. The user will be the only one to assess compliance with the required type of protection and proper connection of the respirator with additional protective equipment.

METHOD OF USING RESPIRATORS:

Before use, check that the selected product is suitable for the existing hazard.

Remove the respirators from the packaging.

Check visually that the respirator is intact (no piercing, no opening, no torn, no rupture, etc.)

Take the respirator from package, put the elastic bands on your ears, covering the mouth and nose. Adjust the nose bridge wire for your nose.

WARNINGS:

*Care should be taken when removing contaminated respirators, so as not to contaminate the user with any hazardous substances. If respirators are contaminated then decontamination procedures should be followed (i.e. decontamination shower) prior to the removal of the respirators. This respirator is not designed for use in extreme environments.

*The determination of suitability of NARKONTEKS products for an application is the final responsibility of the user. This product is recommended for single use application. Upon contamination wear are damage the garment should be removed and appropriately disposed of at the earliest convenience.

*Where NARKONTEKS products are used in conjunction with other PPE, and for full "Type" protection it is necessary to tape cuffs to gloves, ankles to boots, the hood to the respiratory device. The user shall be the sole judge for the correct combination of NARKONTEKS garment and additional PPE.

TRANSPORT, STORAGE AND DISPOSAL:

The product must be transported and stored in the original packaging in dry rooms away from sources of light and heat. If the respirators are not contaminated, they can be disposed of with solid municipal waste. If they have been contaminated, they must be disposed of in accordance with applicable laws and regulations.

USE TIME:

It is suggested to use the product within a period of five years from the date of production written on label Month and year of production: MM/YYYY

EU ATTESTATION OF CONFORMITY

The information on the packaging for the above listed product covers the necessary information stated of Regulation (EU) 2016/425 Module B.


Respiratory protective devices. Filtering half masks to protect against particles. Requirements, testing, marking (EN 149:2001+A1:2009)

CONFORMITY WITH STANDARDS

Protective Disposable Respirator With Valve	FFP2	Result
Total inward leakage	% 11	
Penetration of filter material	% 6 max	
Breathing resistance (inhalation)	30 l/mn 95 l/mn	
	0,7 2,4	
Breathing resistance (exhalation)	160 l/mn	
	3,0	

DIMENSIONS OF PRODUCT

Length	150 +- 5 mm
Width	106 +- 5 mm
Nose Bridge Wire Length	(Min.) 95,0 mm
Elastic Band Length	210 mm +- 5mm



PRODUCT LABEL (EXAMPLE):

Producer: **NARKONTEKS**
TEKSTİL İHR.İTH.SAN. VE TİC.A.Ş.
Dogus Cd, 3/19 Sk, No:12, Begos 35160,
Buca/İZMİR/TURKEY - 0090-232-4408686

Style Code: ALIVE PPE 9020
EN 149:2001+A1:2009



MAINTENANCE AND CLEANING:

